



VASECTOMY KIT UNDER LOCAL ANESTHESIA



o.r. consulting[®]
www.orconsulting.biz

Tél. **+41.22.592.64.74**
4B, chemin de Joran - 1260 Nyon
info@orconsulting.biz

VASECTOMY KIT

VASECTOMY KIT UNDER LOCAL ANESTHESIA

- Insulation Identification and isolation of the deferent channel with the specific ring clamp (1)
- Skin discision incision then dissection of the deferent canal with the fine mosquito pliers (2)
- Externalization and pinching of the deferent canal using Halstead's clamp (3)
- Section of the deferent canal with surgical scissors (4) then coagulation of the light of both ends of the deferent and interposition of tissue between the ends
- The Kits make it possible to perform vasectomies with incisions < 10 mm, local sub-anesthesia, in the operating room as well as in a consultation room or day hospitalization.
- These instruments are made of J1/J2 (Japanese) stainless steel, the hardest and most resistant.
- All our kits are packaged in an autoclavable metal box.
- The usage video is available via the link: <https://youtu.be/ROEMfjqo9VM>



1-A Ring clamp Ø 4mm
for transcutaneous grip of the deferent - 14cm



1-B Ring clamp Ø 6mm
with beveled bits for
transcutaneous grip of the deferent - 14cm



2 Thin and pointed mosquito tweezers
for skin discision (without opening)
and dissection of the deferent canal - 12.5cm



3 Halstead's fine pliers for pinching
and ligation of the deferent - 12.5cm



4 Scissors for section of the deferent - 11cm



5 Dissecting forceps - 14cm

KIT 1

- Ring clamp forceps Ø 4mm, 14cm
- Dissecting sharp forceps, 12.5cm
- Iris scissors, 11cm
- Straight Halstead forceps, 12.5cm

KIT 2

- Ring clamp forceps Ø 4mm, 14cm
- Ring clamp forceps Ø 6mm, 14cm beveled
- Dissecting sharp forceps, 12.5cm
- Iris scissors, 11cm
- Straight Halstead forceps, 12.5cm
- Forceps with jaws, 14cm



o.r. consulting[®]
www.orconsulting.biz
Tél. +41.22.592.64.74
4B, chemin de Joran - 1260 Nyon
info@orconsulting.biz